



Tools for Scientists<sup>®</sup>



# Standard Heating & Controls



**Glas-Col**<sup>®</sup>  
Tools for Scientists<sup>®</sup>

[www.glascol.com](http://www.glascol.com)



# Standard Heating & Controls



Standard "O" Series

## Series "O" is our most popular line of soft heating mantles for round bottom flask

Heating mantles are preferred to Bunsen burners or hot plates for organic liquids. They provide even heating (no hot spots) and maximum surface coverage for better temperature uniformity.

The Series "O" spherical flask mantle covers only the flask's bottom half, letting you see its full contents. It should be used with the Series O supports and a control.

### Fabric Heating Mantles Feature:

- Fabric exterior provides effective heating in a space saving configuration
- Fabric interior to softly nest glass vessels to reduce the chance of thermal shock and damage to glassware
- Withstands 400°C internal operating temperature
- Ideal for bottom holes, viewing ports, and other tubulations
- All 115 volt units are CSA certified
- Electrical: Separable 4 foot 2 wire cord and locking connector. Larger size Series O

mantles have a silicone impregnated exterior for moisture resistance

Catalog Number	Voltage	Flask Capacity (ml)	Watts
100A 0394	115	50	60
100A 0396	115	100	80
100A 0398	115	125	80
100A 0400	115	200	100
100A 0402	115	250	180
100A 0404	115	300	180
100A 0406	115	500	270
100A 0408	115	1,000	380
100A 0410	115	2,000	500
100A 0412	115	3,000	500
100A 0414	115	5,000	600
100A 0416	115	12,000	2@650

Catalog Number	Voltage	Flask Capacity (ml)	Max. Flask Diameter (mm)	Watts
100B TM94	115	50	48	60
100B TM96	115	100	60	80
100B TM98	115	125	70	80
100B TM100	115	200	76	100
100B TM102	115	250	83	180
100B TM104	115	300	86	180
100B TM106	115	500	101	270
100B TM108	115	1,000	137	380
100B TM110	115	2,000	170	500
100B TM112	115	3,000	183	500
100B TM114	115	5,000	220	600
100B TM116	115	12,000	293	2@650
100B TM118	115	22,000	347	2@770
100B TM120	115	50,000	456	3@1,000

## Series "TM" is another popular line of Rigid mantles for round bottom flask


Heating mantles are preferred to Bunsen burners or hot plates for organic liquids. They provide even heating (no hot spots) and maximum surface coverage for better temperature uniformity.

The Series "TM" spherical flask mantle covers only the flask's bottom half, letting you see its full contents. Larger sizes have multiple circuits for ease of temperature control when a flask is less than half full.

### Aluminum Housed Mantles Feature:

- Rigid housing provides strength and durability while supporting the weight of the vessel
- Fabric interior to softly nest glass vessels to reduce the chance of thermal shock and damage to glassware
- Withstands 400°C internal operating temperature
- Feet on smaller sizes designed to provide stability and promote cooler exterior temperatures
- All 115 volt units are CSA certified

mantles have a silicone impregnated exterior for moisture resistance



TM120 shown with terminal box

# Standard Heating & Controls



Catalog Number	Flask Size (ml)	Bottom	Voltage	Watts/Position
100D MX25212	250	Round	120	180
100D MX50212	500	Round	120	270

Replacement Element	Flask Size (ml)	Bottom	Voltage	Watts
100C M102E	250	Round	120	180
100C M106E	500	Round	120	270



## Combo Heating Mantles

Glas-Col has always been the source of mantles for repetitive extracting, refluxing, and distilling procedures in labs of food, textile fiber, water and waste-water, petroleum, and many other industries. Our latest Combo version offers safety and convenience in a multi-place heating mantle.

- Low profile for space-saving convenience
- Easy to replace heating element
- Element containers act as spill containment chambers to capture chemicals from spillage or boil-over
- Clear anodized aluminum cabinet with chemical resistance top
- Offered in two common 3 position configurations: 250ml and 500ml round bottom
- Max Temperature 400°C
- Housing Size: 27"W x 10-3/4"D x 6-1/4"H
- Built in cooling fan to keep housing cool to touch
- Over-temperature safety limit switch
- Lighted power switch for each position

Each Combo is furnished with grid support system and 3 glassware adjustable clamps. Proportional voltage controls are built into each position for constant, steady-state voltage.

**Glassware not included**



# Standard Heating & Controls



### PowerTrol:

This control is small in size and can be used to regulate the input voltage of resistive devices by turning the dial knob. Just plug your resistive load into the receptacle on the back, turn on the switch and adjust the dial. The higher the number on the dial, the hotter your load will get. This unit comes with a 6' attached power cord and grid support bracket, which is ideal for fume hood mounting to maximize bench space.



### Control Features:

- On/Off switch •
- Non-Linear control •  
(scale is for reference only)
- Circuit protection •
- Low profile housing •

Catalog Number	Voltage	Description	Range	Product Weight (lbs)	Watt
104A PL120	120	10A Proportional Voltage Control	from 5-100% of rated voltage	1.5	1200

### DigiTrol II:

This temperature control displays both the setpoint and process temperature and uses the most modern control technology for the best temperature stability. The Auto-tune feature minimizes setpoint overshoot and learns your process. Changing your setpoint is easily done with the interface keys. The controller can also perform ramp rate operations to allow the user to slowly raise the process temperature. This unit comes with a 6' detachable power cord, grid support bracket, which is ideal for fume hood mounting to maximize bench space. This control is a microprocessor-based, digital indicating, automatic temperature control with a single input and a single output. It features an auto-tuning function that allows automatic setting of control parameters with a minimum of user input required. This controls accepts common thermocouples "J", "K", or "T".



The control automatically sets the PID parameters through a "learning" sequence in the auto-tuning mode. PID parameters include proportional band, reset/integral and rate/derivative. User-friendly features include automatic LED indicators to aid in monitoring and setup, as well as dual LED displays for process temperature and set point indication. This control automatically stores all information in a non-volatile memory.

Catalog Number	Voltage	Description	Thermocouple Type	Product Weight (lbs)	Watt
104A PL612	120	PID Control	Universal	3	1800

# Standard Heating & Controls



### RampTrol:

If your heating process requires several steps, this is the control for you. The Ramp/Soak feature of this control allows up to a 40-step profile, which can be configured from the provided software. You can profile Temperature, Time, Hold, Soak and End steps to create the ideal profile for your process. The built in adaptive control technology provides even tighter control for these demanding applications. Several input types are available. This unit comes with 6' detachable power cord, communication port, audible alarm function and grid support bracket, which is ideal for fume hood mounting to maximize bench space.



### RampTrol:

This control is a microprocessor-based, digital indicating, automatic temperature control with a single input and a single output. It features an auto-tuning function that allows automatic setting of control parameters with a minimum of user input required. This control accepts common thermocouples "J", "K", or "T".

The control automatically sets the PID parameters through a "learning" sequence in the auto-tuning mode. PID parameters include proportional band, reset/integral and rate/derivative. User-friendly features include automatic LED indicators to aid in monitoring and setup, as well as dual LED displays for process temperature and set point indication. This control automatically stores all information in a non-volatile memory.

### TwinTrol:

This control is ideal for 2 independent temperature processes. This space saving design combines 2 controls into one package to maximize bench space. This control displays both the setpoint and process temperature and uses the most modern control technology for the best temperature stability. The Auto-tune feature minimizes setpoint overshoot and learns your process. Changing your setpoint is easily done with the interface keys. Several input types are available. This unit comes with an audible alarm function, 6' detachable power cord and grid support bracket, which is ideal for fume hood mounting to maximize bench space. This controls accepts common thermocouples "J", "K", or "T".

The control automatically sets the PID parameters through a "learning" sequence in the auto-tuning mode. PID parameters include proportional band, reset/integral and rate/derivative. User-friendly features include automatic LED indicators to aid in monitoring and setup, as well as dual LED displays for process temperature and set point indication. This control automatically stores all information in a non-volatile memory.

Catalog Number	Voltage	Description	Product Weight (lbs)	Watt
104A PL912	120	PID Control	4	1800



Catalog Number	Voltage	Description	Product Weight (lbs)	Watt
104A TWIN-TROL	120	PID Control (2 Zones)	7	900 per channel

# Standard Heating & Controls



## Specifications & Description

Min Temperature	(°F)-328
Max Temperature	(°F)3214
Min Temperature	(°C)-200
Max Temperature	(°C)1768
Accuracy	0.1% Span
Resolution	0.1°C
Data Logging	Yes
Alarm	Visual/audible by pc
Display	Touchscreen LCD
Input Voltage	100-250 VAC
Thermocouple Types	J, K, T, B, N, S, R & E
Description	12 Channel Scanning Data Logging Benchtop Thermocouple Thermometer



Catalog Number	Voltage	Description	Product Weight (lbs)	Watt
108A TOW-OTP1800	120	PID Control with Limit	7	900 per channel

## Temperature/Limit Control

This control combines a Process control and a Limit control in a space saving design to maximize bench space. The Process control displays both the set-point and process temperature and uses the most modern control technology for the best temperature stability. The Limit control allows the user to input a value that once reached will turn off the power to the load receptacles. The Limit control also displays either a Safe or Fail condition and requires a manual reset. Both the Process and Limit control have audible alarms and separate sensors. A bracket on the back of the control allows for mounting the control on a grid structure. Each load receptacle is rated at 7.5 amps. This control accepts common thermocouples "J", "K", or "T".

## Data Logger

This 12-channel full touchscreen scanning thermometer simultaneously monitors different temperature processes and provides real-time data on your PC with the supplied data acquisition software. Channel configuration can be changed through the software. The thermometer also functions as a data logger, storing readings in internal nonvolatile memory for later download to PC in CSV format.

Visible alarms alert you when a channel logs a temperature exceeding the high/low set point or hysteresis. Thermocouple temperature probes are sold separately.



Catalog Number	Sensor Type	Voltage
108A DTM12	Thermocouple	120/240VAC

# Standard Heating & Controls



# Tools for Scientists<sup>®</sup>



Glas-Col, LLC  
711 Hulman Street, P.O. Box 2128, Terre Haute, IN  
47802-0128 USA  
Phone. (812) 235-6167, Fax. (812) 234-6975  
E-mail: [pinnacle@glascol.com](mailto:pinnacle@glascol.com)  
[www.glascol.com](http://www.glascol.com)



# Glas-Col<sup>®</sup>

Tools for Scientists<sup>®</sup>